

17th INTERNATIONAL FORUM ON MPSoC



# **A Wearable Biomedical Sensing System with Normally-off Computing Architecture**

July 5, 2017


Masahiko Yoshimoto  
Kobe University

# Outline



---

- Background
- Design Issues and Design Features
- Low Power VLSI Design for ECG and Tri-axial Acceleration Monitoring
- Fabrication & Field Testing Results
- Summary

# Population of more than 65 years old



# Suppression of lifestyle disease



Physical activity estimation:  
type,  
intensity,  
duration.


- Exercise
- Intake
- Sleep

Ex.) Diabetes(糖尿病) ↔ Sedentary

Physical activity estimation for sedentary break


# Biomedical data for physical activity estimation

## • Triaxial Acceleration





## • Instantaneous Heart Rate (IHR)

$$= 60 / (\text{latest } RR[s]) \text{ [bpm]}$$



Energy expenditure is estimated by the above two kind of biosensing data.

# Ex) Triaxial acceleration measurement



3-axis  
Acceleration




Activity  
Estimation



Energy  
Expenditure



# Development of wearable bio-medical sensor



- **Requirements**

→ **Low power, Small size and Light weight**

# Design issues

---

- Ultra low power consumption for longer life operation with tiny battery.
- Noise tolerant IHR(Instantaneous Heart Rate) Monitoring at short distance electrode condition to realize wearable small size sensor.



# Design Features

---

- Non-volatile memory and Non-volatile MCU(Micro-Control Unit) for normally-off computing
- Algorithm of noise-tolerant IHR(Instantaneous Heart Rate) extraction

# Normally-off computing

	Frequency component
ECG	0.1 - 150Hz
EEG	0.5 - 60Hz
VEP	0.5 - 60Hz
EMG	few kHz




**Operational clock**

- ✓ CPU
- ✓ MPU


etc...

- **Extremely Low frequency range of Bitral signal**
- **→ Standby power reduction is effective**






**Non-volatile FeRAM was employed to suppress standby leak.**

# VLSI block diagram of bio-medical sensor



# Normally-off computing with Non-volatile memory

## Operating condition





FeRAM has the lowest leak-power feature for bio-medical sensing.

N-words access per second

# 6T-4C non-volatile memory

## 6T-4C shadow memory

6T SRAM + Fe capacitor  
= 6T-4C shadow memory




Active: 6T SRAM

Sleep: Nonvolatile memory

# Design issue of 6T-4C non-volatile memory

Problem

1. Increase of the active energy consumption




2. Increase of the energy consumption during transition period

# Circuit technology for power reduction

---


## 1. Bit-line Non-precharge

→ Reducing the active energy consumption



## 2. Plate-line Charge sharing


→ Reducing the turning-on/off energy consumption




# Bit-line Non-precharge

Proposed method

Connecting a bit-line pair w/ switch



BL: Bit-line  
WL: Word-line  
MC: Memory cell





The energy consumption for charging is reduced



# Plate-line Charge share

Proposed method

When charging PL, reusing an electric charge of previous one




PL: Plate-line


Charging rate is reduced

# Power reduction of 6T4C memory

## 1. Bit-line Non-precharge




## 2. Plate-line Charge sharing




**Improvement of Break-even time**

# VLSI block diagram of bio-medical sensor




# Non-volatile logic

NVFF has capability to retain a logic state w/o power supply.



Ferroelectric capacitor retains data during power off state.

# Basic operation




Logic op: NV-CPU executes instruction in NVRAM at 24MHz.

STORE: Logic state in NVFF is written into FE capacitors.

RECALL: Logic state is recovered from FE capacitors.

# Power dissipation of non-volatile MCU




Non-volatile Flip-Flop circuit (NVFF) is useful to reduce power dissipated in the MCU for vital sensor application.

# Design features







---

- Non-volatile memory and Non-volatile MCU(Micro-Control Unit) for normally-off computing
- Algorithm of noise-tolerant IHR(Instantaneous Heart Rate) extraction

# Various noises in ECG wave



(a) ECG waveform example




(b) Noise problem of threshold approach


(c) Various noises



# QSWT filter for noise suppression




(a) Block diagram




(b) Frequency characteristics


# Algorithm of noise torelant IHR extraction



(a) Coarse QRS search





(b) Fine QRS search



(c) Template matching


Coarse-fine QRS template generation and template matching with QRS prediction.

# Success rate evaluation of IHR extractor




The IHR extractor can also suppress motion artifact and muscle artifacts.

# Photo of VLSI wearable bio-medical sensor




Technology	0.13 $\mu\text{m}$ CMOS	
Supply voltage	1.2 V (AFE, ADC, Logic, Memory)	
	3.0 V (32.768-kHz OSC and I/O)	
Chip area	3.7 $\times$ 4.3 mm <sup>2</sup>	
Frequency	24 MHz (for processor)	
	32 kHz (for other blocks)	
Processor	32-bit Cortex M0 (with NVFF)	
On chip memory	16-KByte 6T-4C NVRAM	
AFE	Gain	54 dB
	Bandwidth	700 Hz
	CMRR	73 dB
ADC	Resolution	8 bit
	Current	0.23 $\mu\text{A}$ @128 S/s, 1.0 $\mu\text{A}$ @1 kS/s


# Experimental wearable bio-medical sensor



# Measurement results




(a) Total power




(b) IHR extractor comparison

# Power reduction in sensor module



# Technical position of the bio-medical sensor





# Summary

---

- The wearable bio-medical sensor has been developed, which features heart rate monitoring and tri-axial acceleration using newly developed low power SoC.
- Non-volatile MCU for normally-off computing and noise-tolerant IHR detection algorithm have been employed for ECG-SoC design.
- The fully integrated ECG-SoC consumes  $6.14\mu\text{A}$  for ECG monitoring and the sensor system dissipates  $20\mu\text{A}$ , allowing 2-weeks continuous sensing only using a 10-mAh thin-type lithium-ion battery.
- The activity classification using data from field testing attains ~90% accuracy.

# Acknowledgement

---

Author would like to thank  
OMRON HEALTHCARE Co., Ltd for  
algorithm and software development and  
ROHM SEMICONDUCTOR for fabrication of  
SoC and experimental module.

This research was supported by Ministry of  
Economy, Trade and Industry (METI) and the  
New Energy and Industrial Technology  
Development Organization (NEDO).